

# CEA antibody

Catalog Number: 109242

## Product name

CEA antibody

## Specificity

Human; other species not tested.

## Antibody description

CEA Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue, human colon tissue, human lung cancer tissue, human pancreas cancer tissue, human stomach cancer tissue.

## Preparation

This antibody was obtained by immunization of CEA recombinant protein (Accession Number: NM\_004363). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:109242(CEA antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:109242(CEA antibody) at dilution of 1:50 (under 40x lens)