

# CEACAM3-Specific antibody



Catalog Number: 109244

## Product name

CEACAM3-Specific antibody

## Specificity

Human; other species not tested.

## Antibody description

CEACAM3-Specific Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue, human stomach cancer tissue.

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_001815). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:109244(CEACAM3-Specific antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:109244(CEACAM3-Specific antibody) at dilution of 1:50 (under 40x lens)