

CENPF antibody

Catalog Number: 109252

Product name

CENPF antibody

Specificity

Human; other species not tested.

Antibody description

CENPF Rabbit Polyclonal antibody. Positive IHC detected in human kidney tissue, human breast cancer tissue. Positive IF detected in HeLa cells.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_016343). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, IF


Dilutions

Recommended Dilution:


IHC: 1:20-1:200

IF: 1:20-1:200


Validations



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:109252(CENPF Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:109252(CENPF Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HeLa cells using Catalog No:109252(CENPF Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)