

Centrin 3 antibody



Catalog Number: 109261

Product name

Centrin 3 antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat; other species not tested.

Applications

ELISA, WB

Antibody description

Centrin 3 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, mouse ovary tissue. Observed molecular weight by Western-blot: 23 kDa

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of Centrin 3 recombinant protein (Accession Number: NM_004365). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:109261(CETN3 antibody) at dilution of 1:300

Clonality