

CEP110 antibody

Catalog Number: 109263

Product name

CEP110 antibody

Immunogen

[Human CEP110 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

CEP110 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

Preparation

This antibody was obtained by immunization of CEP110 recombinant protein (Accession Number: NM_007018). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:109263(CEP110 Antibody) at dilution of 1:700