

CEP192 antibody



Catalog Number: 109269

Product name

CEP192 antibody

Specificity

Human; other species not tested.

Antibody description

CEP192 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 250-300 kDA

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: BC144481). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, IF


Dilutions

Recommended Dilution:


IP: 1:200-1:1000

IF: 1:20-1:200

Validations



IP result of anti-CEP192 (Catalog No:109269 for IP and Detection) with HeLa cells.



Immunofluorescent analysis of HeLa cells using Catalog No:109269(CEP192 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)