

CHRAC1 antibody

Catalog Number: 109275

Product name

CHRAC1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CHRAC1 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, L02 cells. Positive IHC detected in human pancreas cancer tissue. Observed molecular weight by Western-blot: 15 kDa

Preparation

This antibody was obtained by immunization of CHRAC1 recombinant protein (Accession Number: NM_017444). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:109275(CHRAC1 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:109275(CHRAC1 antibody) at dilution of 1:50 (under 10x lens)