

CHRNA6 antibody



Catalog Number: 109281

Product name

CHRNA6 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CHRNA6 Rabbit Polyclonal antibody. Positive WB detected in mouse thymus tissue. Observed molecular weight by Western-blot: 57 kDa

Preparation

This antibody was obtained by immunization of CHRNA6 recombinant protein (Accession Number: XM_047422396). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:109281(CHRNA6 antibody) at dilution of 1:500