

CHRNB4 antibody



Catalog Number: 109286

Product name

CHRNB4 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

CHRNB4 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, HeLa cells, mouse testis tissue. Positive FC detected in SH-SY5Y cells. Observed molecular weight by Western-blot: 48-60 kDa

Preparation

This antibody was obtained by immunization of CHRNB4 recombinant protein (Accession Number: NM_000750). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, FC

Dilutions


Recommended Dilution:

WB: 1:1000-1:10000

Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:109286(CHRNB4 antibody) at dilution of



1X10⁶ SH-SY5Y cells were stained with 0.2ug CHRNB4 antibody (Catalog No:109286, red) and control antibody (blue). Fixed with 4% PFA blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.