

CHST15 antibody



Catalog Number: 109296

Product name

CHST15 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CHST15 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue, Jurkat cells, mouse spleen tissue, Raji cells. Positive IP detected in Raji cells. Observed molecular weight by Western-blot: 49 kDa, 69 kDa

Preparation

This antibody was obtained by immunization of CHST15 recombinant protein (Accession Number: BC050540). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:109296(CHST15 antibody) at dilution of 1:500



IP Result of anti-CHST15 (IP:Catalog No:109296, 4ug; Detection:Catalog No:109296 1:500) with Raji cells lysate 3200ug.