

# CHST15 antibody

Catalog Number: 109296

## Product name

CHST15 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CHST15 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue, Jurkat cells, mouse spleen tissue, Raji cells. Positive IP detected in Raji cells. Observed molecular weight by Western-blot: 49 kDa, 69 kDa

## Preparation

This antibody was obtained by immunization of CHST15 recombinant protein (Accession Number: BC050540). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:109296(CHST15 antibody) at dilution of 1:500



IP Result of anti-CHST15 (IP:Catalog No:109296, 4ug; Detection:Catalog No:109296 1:500) with Raji cells lysate 3200ug.