

# CHST2 antibody



Catalog Number: 109297

## Product name

CHST2 antibody

Polyclonal

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse; other species not tested.

## Applications

ELISA, WB

## Antibody description

CHST2 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 50-52, 58-60 kDa

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

## Preparation

This antibody was obtained by immunization of CHST2 recombinant protein (Accession Number: NM\_004267). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109297(CHST2 Antibody) at dilution of 1:600

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality