

cIAP1 antibody

Catalog Number: 109308

Product name

cIAP1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

cIAP1 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, Jurkat cells, mouse brain tissue, mouse liver tissue, mouse testis tissue. Positive IHC detected in human testis tissue, human brain tissue, human kidney tissue, human ovary tissue, human pancreas tissue, human placenta tissue, human spleen tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 70 kDa, 55-60 kDa

Preparation

This antibody was obtained by immunization of cIAP1 recombinant protein (Accession Number: NM_001166). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF, WB, IHC

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


IHC: 1:20-1:200

IF: 1:10-1:100

Validations




mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:109308(cIAP1 antibody) at dilution of 1:1000




Immunohistochemical of paraffin-embedded human testis using Catalog No:109308(cIAP1 antibody) at dilution of 1:50 (under 10x lens)

cIAP1 antibody

Catalog Number: 109308



Immunohistochemical of paraffin-embedded human testis using Catalog No:109308(cIAP1 antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using BIRC2 antibody Catalog No:109308 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).