

Product name

CIDEA antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CIDEA Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue, mouse brain tissue. Observed molecular weight by Western-blot: 25kd

Preparation

This antibody was obtained by immunization of CIDEA recombinant protein (Accession Number: NM_001318383). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations

mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:109313(CIDEA antibody) at dilution of 1:100