

CIDEC antibody



Catalog Number: 109314

Product name

CIDEC antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CIDEC Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells, HeLa cells, L02 cells. Observed molecular weight by Western-blot: 29 kDa; 48 kDa

Preparation

This antibody was obtained by immunization of CIDEC recombinant protein (Accession Number: XM_047448685). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:109314(CIDEC Antibody) at dilution of 1:1000