

CKMT2 antibody



Catalog Number: 109331

Product name

CKMT2 antibody

IHC: 1:20-1:200


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

CKMT2 Rabbit Polyclonal antibody. Positive WB detected in human skeletal muscle tissue, human heart tissue, Jurkat cells. Positive IHC detected in human skeletal muscle tissue, human heart tissue. Observed molecular weight by Western-blot: 48 kDa




Preparation

This antibody was obtained by immunization of CKMT2 recombinant protein (Accession Number: NM_001099736). Purification method: Antigen affinity purified.

human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:109331(CKMT2 antibody) at dilution of 1:800

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal


Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using Catalog No:109331(CKMT2 Antibody) at dilution of 1:50 (under 10x lens)

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC



Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using Catalog No:109331(CKMT2 Antibody) at dilution of 1:50

CKMT2 antibody



Catalog Number: 109331

(under 40x lens)