

Claudin 22 antibody



Catalog Number: 109339

Product name

Claudin 22 antibody

Specificity

Human; other species not tested.

Antibody description

Claudin 22 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, A431 cells. Observed molecular weight by Western-blot: 25-30 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001111319). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109339(CLDN22 antibody) at dilution of 1:500