

CLEC10A/CD301 antibody



Catalog Number: 109346

Product name

Polyclonal

CLEC10A/CD301 antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

CLEC10A/CD301 Rabbit Polyclonal antibody.
Positive WB detected in BxPC-3 cells. Observed
molecular weight by Western-blot: 32 kDa, 28 kDa

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of
CLEC10A/CD301 recombinant protein (Accession
Number: NM_182906). Purification method:
Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



BxPC-3 cells were subjected to SDS PAGE
followed by western blot with Catalog
No:109346(CLEC10A antibody) at dilution of
1:400