

CLEC14A antibody



Catalog Number: 109348

Product name

CLEC14A antibody

Immunogen

[Human CLEC14A Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

CLEC14A Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 52 kDa

Preparation

This antibody was obtained by immunization of CLEC14A recombinant protein (Accession Number: NM_175060). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:109348(CLEC14A antibody) at dilution of 1:1000