

# CLEC16A antibody



Catalog Number: 109349

## Product name

CLEC16A antibody

## Immunogen

[Human CLEC16A Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

CLEC16A Rabbit Polyclonal antibody. Positive WB detected in mouse spleen tissue. Observed molecular weight by Western-blot: 104-110 kDa

## Preparation

This antibody was obtained by immunization of CLEC16A recombinant protein (Accession Number: NM\_001243403). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:109349(CLEC16A Antibody) at dilution of 1:600