

CNTN4 antibody

Catalog Number: 109363

Product name

CNTN4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CNTN4 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 113kd

Preparation

This antibody was obtained by immunization of CNTN4 recombinant protein (Accession Number: NM_001206956). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109363(CNTN4 antibody) at dilution of 1:300