

COG1 antibody

Catalog Number: 109367

Product name

COG1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

COG1 Rabbit Polyclonal antibody. Positive IHC detected in human ovary tumor tissue. Positive WB detected in mouse testis tissue. Observed molecular weight by Western-blot: 150kd+100kd+65kd

Preparation

This antibody was obtained by immunization of COG1 recombinant protein (Accession Number: NM_018714). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:109367(COG1 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:109367(COG1 antibody) at dilution of 1:100 (under 10x lens)