

# COG6 antibody



Catalog Number: 109369

## Product name

COG6 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

COG6 Rabbit Polyclonal antibody. Positive FC detected in MCF-7 cells. Positive IF detected in MCF-7 cells, Hela cells. Positive WB detected in human brain tissue, NIH/3T3 cells. Observed molecular weight by Western-blot: 73kd

## Preparation

This antibody was obtained by immunization of COG6 recombinant protein (Accession Number: NM\_020751). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF, FC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100


## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109369(COG6 antibody) at dilution of 1:300



Immunofluorescent analysis of MCF-7 cells, using COG6 antibody Catalog No:109369 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).



1X10<sup>6</sup> MCF-7 cells were stained with 0.2ug

# COG6 antibody



Catalog Number: 109369

---

COG6 antibody (Catalog No:109369, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-

conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.