

COLEC11 antibody



Catalog Number: 109374

Product name

COLEC11 antibody

Specificity

Human, Mouse, Rat, Zebrafish; other species not tested.

Antibody description

COLEC11 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, human blood tissue, rat liver tissue. Positive IP detected in mouse liver tissue. Positive IF detected in ATCD5 cells. Observed molecular weight by Western-blot: 34 kDa

Preparation

This antibody was obtained by immunization of COLEC11 recombinant protein (Accession Number: NM_024027). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


IP: 1:200-1:2000

IF: 1:20-1:200

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:109374(COLEC11 antibody) at dilution of 1:500



IF result(confocal image) of CL-K1 in ATCD5 cell show Golgi localization from Prof. Philip L. Beales.



IP Result of anti-COLEC11 (IP:Catalog No:109374, 4ug; Detection:Catalog No:109374 1:500) with mouse liver tissue lysate 6400ug.