

# Collagen Type II antibody



Catalog Number: 109376

## Product name

Collagen Type II antibody

## Specificity

Human; other species not tested.

## Antibody description

Collagen Type II Rabbit Polyclonal antibody.  
Positive WB detected in transfected cells. Positive  
IHC detected in human lung tissue.

## Preparation

This antibody was obtained by immunization of  
Collagen Type II recombinant protein (Accession  
Number: BC007252). Purification method: Antigen  
affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH  
7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:250-1:1000

IHC: 1:20-1:200

## Validations




WB result of anti-COL2A1 (Catalog No:109376)  
with overexpressed protein COL2A1(1-268aa)-



GFP.

Immunohistochemistry of paraffin-embedded  
human lung tissue slide using Catalog  
No:109376(COL2A1 Antibody) at dilution of 1:50  
(under 10x lens)



Immunohistochemistry of paraffin-embedded  
human lung tissue slide using Catalog  
No:109376(COL2A1 Antibody) at dilution of 1:50  
(under 40x lens)