

Collagen Type IV antibody



Catalog Number: 109379

Product name

Collagen Type IV antibody

Specificity

Human, Rat; other species not tested.

Antibody description

Collagen Type IV Rabbit Polyclonal antibody. Positive IHC detected in human heart tissue, human kidney tissue, rat brain tissue.

Preparation

This antibody was obtained by immunization of Collagen Type IV recombinant protein (Accession Number: NM_001845). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human heart slide using Catalog No:109379(COL4A1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human heart slide using Catalog No:109379(COL4A1 Antibody) at dilution of 1:50 (under 40x lens)