

CLIC4 antibody



Catalog Number: 109390

Product name

CLIC4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CLIC4 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, human kidney tissue, mouse ovary tissue. Observed molecular weight by Western-blot: 29 kDa

Preparation

This antibody was obtained by immunization of CLIC4 recombinant protein (Accession Number: NM_013943). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:109390(CLIC4 Antibody) at dilution of 1:600