

# CLIP-115/CYLN2 antibody



Catalog Number: 109393

## Product name

CLIP-115/CYLN2 antibody

## Immunogen

[Human CLIP2 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CLIP-115/CYLN2 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, rat brain tissue. Observed molecular weight by Western-blot: 110-120 kDa

## Preparation

This antibody was obtained by immunization of CLIP-115/CYLN2 recombinant protein (Accession Number: NM\_032421). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:1000-1:10000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109393(CLIP2 Antibody) at dilution of 1:3000