

# CLIP4 antibody



Catalog Number: 109394

## Product name

CLIP4 antibody

## Specificity

Human; other species not tested.

## Antibody description

CLIP4 Rabbit Polyclonal antibody. Positive IF detected in HUVEC cells. Positive WB detected in HUVEC cells, Y79 cells. Observed molecular weight by Western-blot: 65-75 kDa

## Preparation

This antibody was obtained by immunization of CLIP4 recombinant protein (Accession Number: BC015310). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:10-1:100


## Validations



HUVEC cells were subjected to SDS PAGE followed by western blot with Catalog No:109394(CLIP4 Antibody) at dilution of 1:300



Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:109394(CLIP4 Antibody) at dilution of 1:600



Immunofluorescent analysis of (-20oc Ethanol) fixed HUVEC cells using Catalog No:109394(CLIP4 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)