

CLRN2 antibody



Catalog Number: 109403

Product name

CLRN2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

CLRN2 Rabbit Polyclonal antibody. Positive WB detected in mouse eye tissue. Observed molecular weight by Western-blot: 25 kDa

Preparation

This antibody was obtained by immunization of CLRN2 recombinant protein (Accession Number: NM_001079827). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse eye tissue were subjected to SDS PAGE followed by western blot with Catalog No:109403(CLRN2 Antibody) at dilution of 1:300