

# c-Myb antibody

Catalog Number: 109417

## Product name

c-Myb antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

c-Myb Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells, Jurkat cells. Positive IP detected in HL-60 cells. Observed molecular weight by Western-blot: 72-80kd

## Preparation

This antibody was obtained by immunization of c-Myb recombinant protein (Accession Number: NM\_005375). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

## Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:109417(MYB antibody) at dilution of 1:400



IP Result of anti-MYB (IP:Catalog No:109417, 4ug; Detection:Catalog No:109417 1:300) with HL-60 cells lysate 4000ug.