

CNOT10 antibody

Catalog Number: 109432

Product name

CNOT10 antibody

WB: 1:500-1:5000

IP: 1:200-1:1000

Specificity

Human; other species not tested.


IHC: 1:20-1:200

IF: 1:10-1:100

Antibody description

CNOT10 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in human brain tissue, human testis tissue. Positive IHC detected in human colon cancer tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 82 kDa

Validations



Preparation

This antibody was obtained by immunization of CNOT10 recombinant protein (Accession Number: NM_001256741). Purification method: Antigen affinity purified.

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109432(CNOT10 antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

ELISA, WB, IHC, IF, IP

Dilutions


Recommended Dilution:




Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:109432(CNOT10 antibody) at dilution of 1:50 (under 10x lens)

CNOT10 antibody


Catalog Number: 109432



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:109432(CNOT10 antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using CNOT10 antibody Catalog No:109432 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

IP Result of anti-CNOT10 (IP:Catalog No:109432, 4ug; Detection:Catalog No:109432 1:300) with mouse brain tissue lysate 5200ug.