

Collagen Type XVI antibody



Catalog Number: 109441

Product name

Collagen Type XVI antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Collagen Type XVI Rabbit Polyclonal antibody.
Positive WB detected in human placenta tissue.
Observed molecular weight by Western-blot:
158kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001856).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:109441(COL16A1-Specific antibody) at dilution of 1:1500