

Collagen Type XXIII antibody



Catalog Number: 109442

Product name

Collagen Type XXIII antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Collagen Type XXIII Rabbit Polyclonal antibody.
Positive WB detected in mouse thymus tissue.
Observed molecular weight by Western-blot:
68-75kd

Preparation

This antibody was obtained by immunization of Collagen Type XXIII recombinant protein (Accession Number: BC042428). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:109442(COL23A1 antibody) at dilution of 1:800