

Complement factor H antibody



Catalog Number: 109447

Product name

Complement factor H antibody

Specificity

Human; other species not tested.

Antibody description

Complement factor H Rabbit Polyclonal antibody.
Positive WB detected in human blood tissue.
Positive IP detected in human plasma tissue.
Observed molecular weight by Western-blot: 150 kDa

Preparation

This antibody was obtained by immunization of Complement factor H recombinant protein (Accession Number: NM_001014975). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

Validations



human blood tissue were subjected to SDS PAGE followed by western blot with Catalog No:109447(CFH Antibody) at dilution of 1:1000



IP Result of anti-CFH (IP:Catalog No:109447, 4ug; Detection:Catalog No:109447 1:1000) with human plasma tissue lysate 4000ug.