

# Connexin-46 antibody



Catalog Number: 109454

## Product name

Connexin-46 antibody

## Specificity

Human; other species not tested.

## Antibody description

Connexin-46 Rabbit Polyclonal antibody. Positive WB detected in A431 cells. Observed molecular weight by Western-blot: 46 kDa

## Preparation

This antibody was obtained by immunization of Connexin-46 recombinant protein (Accession Number: XM\_011535048). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:109454(GJA3 Antibody) at dilution of 1:600