

Connexin-46 antibody



Catalog Number: 109454

Product name

Connexin-46 antibody

Specificity

Human; other species not tested.

Antibody description

Connexin-46 Rabbit Polyclonal antibody. Positive WB detected in A431 cells. Observed molecular weight by Western-blot: 46 kDa

Preparation

This antibody was obtained by immunization of Connexin-46 recombinant protein (Accession Number: XM_011535048). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:109454(GJA3 Antibody) at dilution of 1:600