

COX10 antibody



Catalog Number: 109482

Product name

COX10 antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat; other species not tested.

Applications

ELISA, WB

Antibody description

COX10 Rabbit Polyclonal antibody. Positive WB detected in A549 cells, mouse colon tissue, rat colon tissue. Observed molecular weight by Western-blot: 51 kDa

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Preparation

This antibody was obtained by immunization of COX10 recombinant protein (Accession Number: NM_001303). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:109482(COX10 antibody) at dilution of 1:300

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality