

# COX11 antibody

Catalog Number: 109483

## Product name

COX11 antibody

IP: 1:500-1:5000

IHC: 1:20-1:200

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

COX11 Rabbit Polyclonal antibody. Positive IP detected in mouse liver tissue. Positive WB detected in mouse liver tissue, mouse brain tissue. Positive IHC detected in human pancreas cancer tissue. Observed molecular weight by Western-blot: 31 kDa

## Preparation

This antibody was obtained by immunization of COX11 recombinant protein (Accession Number: XM\_017024193). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IHC, IP

## Dilutions


Recommended Dilution:

WB: 1:1000-1:10000


## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:109483(COX11 antibody) at dilution of 1:3000



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:109483(COX11 antibody) at dilution of 1:50



(under 10x lens)

IP Result of anti-COX11 (IP:Catalog No:109483,

# COX11 antibody



Catalog Number: 109483

---

4ug; Detection:Catalog No:109483 1:1000) with mouse liver tissue lysate 4400ug.