

COX16 antibody



Catalog Number: 109485

Product name

COX16 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

COX16 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, A375 cells, HepG2 cells, human heart tissue. Observed molecular weight by Western-blot: 12kd

Preparation

This antibody was obtained by immunization of COX16 recombinant protein (Accession Number: NM_016468). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:109485(COX16 antibody) at dilution of 1:1000