

# COX17 antibody

Catalog Number: 109486

## Product name

COX17 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

COX17 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, rat heart tissue. Positive IHC detected in human lymphoma tissue. Observed molecular weight by Western-blot: 7kd

## Preparation

This antibody was obtained by immunization of COX17 recombinant protein (Accession Number: NM\_005694). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:109486(COX17 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:109486(COX17 antibody) at dilution of 1:50 (under 10x lens)