

COX18 antibody



Catalog Number: 109487

Product name

COX18 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

COX18 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in mouse liver tissue, HeLa cells, NIH/3T3 cells. Observed molecular weight by Western-blot: 37kd-40,28-30kd

Preparation

This antibody was obtained by immunization of COX18 recombinant protein (Accession Number: NM_001033760). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



Immunofluorescent analysis of HepG2 cells, using COX18 antibody Catalog No:109487 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:109487(COX18 Antibody) at dilution of 1:600