

COX6B1 antibody



Catalog Number: 109494

Product name

COX6B1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

COX6B1 Rabbit Polyclonal antibody. Positive WB detected in human liver tissue, human heart tissue. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 10kd

Preparation

This antibody was obtained by immunization of COX6B1 recombinant protein (Accession Number: NM_001863). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:109494(COX6B1 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:109494(COX6B1 antibody) at dilution of 1:50 (under 10x lens)