

# COX7B antibody



Catalog Number: 109500

## Product name

COX7B antibody

## Specificity

Human; other species not tested.

## Antibody description

COX7B Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue. Positive WB detected in human skeletal muscle tissue, human brain tissue, human liver tissue. Observed molecular weight by Western-blot: 9 kDa

## Preparation

This antibody was obtained by immunization of COX7B recombinant protein (Accession Number: NM\_001866). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

## Validations



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:109500(COX7B antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:109500(COX7B antibody) at dilution of 1:100 (under 10x lens)