

COX7B antibody



Catalog Number: 109500

Product name

COX7B antibody

Specificity

Human; other species not tested.

Antibody description

COX7B Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue. Positive WB detected in human skeletal muscle tissue, human brain tissue, human liver tissue. Observed molecular weight by Western-blot: 9 kDa

Preparation

This antibody was obtained by immunization of COX7B recombinant protein (Accession Number: NM_001866). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:109500(COX7B antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:109500(COX7B antibody) at dilution of 1:100 (under 10x lens)