

CPEB1 antibody

Catalog Number: 109506

Product name

CPEB1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CPEB1 Rabbit Polyclonal antibody. Positive IHC detected in human liver cancer tissue, human breast cancer tissue, human colon cancer tissue. Positive WB detected in A375 cells, human heart tissue. Observed molecular weight by Western-blot: 54kd

Preparation

This antibody was obtained by immunization of CPEB1 recombinant protein (Accession Number: NM_001387073). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:109506(CPEB1 antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:109506(CPEB1 Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:109506(CPEB1 Antibody) at dilution of 1:50 (under 40x lens)

CPEB1 antibody



Catalog Number: 109506

No:109506(CPEB1 Antibody) at dilution of 1:50 (under 40x lens)