

# CREB3L2 antibody



Catalog Number: 109548

## Product name

CREB3L2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CREB3L2 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, mouse liver tissue, mouse testis tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 60-70 kDa

## Preparation

This antibody was obtained by immunization of CREB3L2 recombinant protein (Accession Number: NM\_001253775). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:109548(CREB3L2 Antibody) at dilution of 1:600



Immunofluorescent analysis of (-20oc Ethanol) fixed HepG2 cells using Catalog No:109548(CREB3L2 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)