

CREB3L4 antibody



Catalog Number: 109549

Product name

CREB3L4 antibody

Specificity

human; other species not tested.

Antibody description

CREB3L4 Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells. Observed molecular weight by Western-blot: 50kd

Preparation

This antibody was obtained by immunization of CREB3L4 recombinant protein (Accession Number: NM_001255978). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:109549(CREB3L4 antibody) at dilution of 1:500