

CRHR2 antibody



Catalog Number: 109558

Product name

Polyclonal

CRHR2 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

CRHR2 Rabbit Polyclonal antibody. Positive WB detected in rat skeletal muscle tissue, mouse heart tissue. Observed molecular weight by Western-blot: 48 kDa

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001883). Purification method: Antigen affinity purified.

Formulation


PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Validations



rat skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:109558(CRHR2 Antibody) at dilution of 1:600