

CRK antibody



Catalog Number: 109565

Product name

CRK antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

Human, Mouse; other species not tested.

Applications

ELISA, WB

Antibody description

CRK Rabbit Polyclonal antibody. Positive WB detected in A549 cells, NIH/3T3 cells, PC-3 cells. Observed molecular weight by Western-blot: 34 kDa

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of CRK recombinant protein (Accession Number: NM_016823). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:109565(CRK Antibody) at dilution of 1:1000

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality